	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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Maxpro™303 medium is a chemical defined basal medium with no animal-derived ingredients, and no protein ingredients, which is suitable for batch culture, fed-batch culture and perfusion culture in the development and production of therapeutic protein products by Chinese hamster ovary (CHO). Maxpro™303 medium contains no L-glutamine. It is suitable for the culture of different cell lines, such as CHOK1, CHOK1SV and CHOS, using GS screening system.

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Maxpro™303 medium is specially designed to improve the growth and production performance of CHO cells. After cell recovery, it is recommended to use Maxgrow™201 or customer customized culture medium as the seed culture medium. Seed cells at  $0.3\sim 0.8\times 10^6$  cells/mL density for passage amplification; In the N-1st generation, continue amplification using Maxpro™303 as the basal medium, and amplify to the required seed volume; In the Nth generation, Maxpro™303 was used as the basal medium and the inoculation density was controlled between  $0.5\sim 2.0\times 10^6$  cells/mL for fed batch culture. It is recommended to use it with Maxfeed™403A1 feed medium.

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This product does not contain insulin. If insulin is added, the recommended concentration is  $0.5\sim 10\text{mg/L}$ .

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Suitable for powder packaging (take 1L as an example)

1. 准备配制体积 80%左右的超纯水 (20~30℃) ;  
Prepare 80% volume ultrapure water (20 ~ 30℃);
2. 加入 Maxpro™ 303 培养基粉末 21.23g, 搅拌 30 min;  
Add 21.23 g of Maxpro™ 303 medium powder and stir for 30 min;
3. 加入碳酸氢钠 1.49g;  
Add 1.49g sodium bicarbonate;
4. 调节 pH 至 7.00±0.1;  
Adjust pH to 7.00±0.1;
5. 定容, 搅拌 5~10 min;  
Constant volume, stir for 5~10 min;
6. 采用 NaCl 调节渗透压至 280~300 mOsm/kg;  
Use NaCl to adjust the osmotic pressure to 280-300 mOsm/kg;  
备注: 1g 氯化钠约等于 33 mOsm/kg  
Note: 1g of NaCl is approximately equal to 33 mOsm/kg
7. 用 0.22µm 过滤器除菌过滤。  
Use 0.22µm filter to sterilization and filtration.

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Shanghai production base, powder packaging: 2℃ to 8℃, protect from light; validity period: 24 months.

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Wuxi production base, powder packaging: 2℃ to 8℃, protect from light; retest date: 24 months.

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地址 Address: 上海市松江区民强路 1525 弄 30 幢 Building 30, Lane 1525, Minqiang Road, Songjiang District, Shanghai

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